

SHENTEK

CAR/TCR Transgene Copy Number Quantitation Kit

User Guide

This is a watermark for the trial version, register to get the full one!

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.

Remove Watermark Now

Version: A/1

For Research Use Only

Product No.: SK030221CA100

Reagents for 100 Reactions

Huzhou Shenke Biotechnology Co., Ltd

(IMPORTANT: Please read this document carefully before experiment.)

1. Product information

■ Product description

SHENTEK® CAR/TCR Transgene Copy Number Quantitation Kit is designed to quantify the copy number of target genes in human cell products derived from common HIV-1 lentiviral vectors, such as CAR or TCR genes integrated in CAR-T and TCR-T cell genomes.

This kit utilizes multiplex qPCR technique to detect DNA sequences related to the integration or expression function on the transfer plasmid and Single Copy Gene (SCG) in human cells. Average copies of target gene per cell can be calculated in the samples, such as CAR or TCR gene copy numbers. This kit is rapid, specific, reliable, and universal for lentiviral vector applications. The reference standard in

This is a watermark for the trial version, register to get the full one!

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.

Remove Watermark Now

Table 1. Kit components and storage

Reagent	Part No.	Quantity	Storage
qPCR Reaction Buffer	NNB001	850 μ L \times 2 tubes	-20°C, protect from light
CAR/TCR Primer&Probe MIX	NNC032	300 μ L \times 1 tube	
IPC MIX	NNC066	150 μ L \times 1 tube	
DNA Dilution Buffer (DDB)	NND001	1.5 mL \times 3 tubes	-20°C
CAR/TCR Control	NNA022	50 μ L \times 1 tube	

The kit components can be stored at appropriate conditions for up to 24 months.

Please check the expiration date on the labels.

■ Applied instruments, including but not limited to the following

- SHENTEK-96S Real-Time PCR System

- 7500 Real-Time PCR System
- LineGene 9600plus Real-Time PCR System
- CFX96 Real-Time PCR System

■ Required materials not included in the kit

- Nonstick, RNase-free, Low Retention Microfuge Tubes, 1.5mL
- Nonstick, Low Retention Tips: 1000 μ L, 100 μ L and 10 μ L
- 96-well qPCR plates with sealing film or PCR 8-strip tubes with caps

■ Related equipment

- Benchtop microcentrifuge
- Vortex mixer

➤ Microplate shaker

➤ Micropipettes 1000 μ L, 100 μ L and 10 μ L

➤ Real-time PCR System

This is a watermark for the trial version, register to get the full one!

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.

[Remove Watermark Now](#)

Sample preparation



qPCR reaction mix preparation



qPCR amplification



Data analysis

2. Methods

■ Experiment preparation

1. Wear appropriate protective eyewear, mask, clothing and gloves.
2. Irradiate the tabletop, micropipettes and tubes with UV for 30 minutes, and disinfect with 75% ethanol.

3. Thaw the kit completely at 2-8°C or melt on ice.

■ Dilution of CAR/TCR Control and preparation of standard curve

Please check the concentration on the label of the CAR/TCR Control tube prior to dilution. Conversion of copy numbers according to the formula:

$$\text{Plasmid copy numbers (copies/}\mu\text{L)} = 6.02 \times 10^{14} \times \text{Plasmid concentration (ng/}\mu\text{L)} / (\text{Number of plasmid bases} \times 660)$$

The copy number of CAR/TCR and SCG is 2.33×10^8 copies/ μL

Prepare CAR/TCR Control solution with DNA Dilution Buffer (DDB) following the serial dilution procedure below:

1. Thaw CAR/TCR Control and DDB completely at 2-8°C or melt on ice. Vortex to mix well and quickly spin down the reagents for 3-5 seconds in a microcentrifuge, and repeat for 3 times.

This is a watermark for the trial version, register to get the full one!

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.

Remove Watermark Now

4. Add 90 μL DDB to each tube of ST1, ST2, ST3, ST4 and ST5.

5. Perform the serial dilutions according to Table 2:

Table 2. Dilution for CAR/TCR Control

Serial dilution tube	Dilution	Conc. (copies/ μL)	
		CAR/TCR	SCG
ST1	10 μL ST0 + 90 μL DDB	2.33×10^6	2.33×10^6
ST2	10 μL ST1 + 90 μL DDB	2.33×10^5	2.33×10^5
ST3	10 μL ST2 + 90 μL DDB	2.33×10^4	2.33×10^4
ST4	10 μL ST3 + 90 μL DDB	2.33×10^3	2.33×10^3
ST5	10 μL ST4 + 90 μL DDB	2.33×10^2	2.33×10^2

- The remaining unused DDB need to be stored at 2-8°C. If the solution is cloudy or contains precipitates, heat at 37°C until it clear.

■ Sample preparation

➤ Negative Control Sample (NCS) Preparation

Add 100 µL of DDB to a new 1.5 mL microcentrifuge tube, and label as NCS.

NCS and samples should be prepared in same way for DNA extraction.

■ qPCR MIX preparation

1. Determine the number of reaction wells based on the standard curve, with the number of test samples and control samples. Generally, triplicates are tested for each sample.

Number of reaction wells = (5 standard points on the standard curve + 1 NTC + 1 NCS + test samples) × 3

2. Prepare qPCR MIX according to the number of reaction wells.

Table 3. qPCR MIX preparation

This is a watermark for the trial version, register to get the full one!

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.

Remove Watermark Now

3. Mix thoroughly and place on ice, aliquot 20 µL/well into 96-well qPCR plate or PCR 8-strip tubes.

■ qPCR Reaction MIX preparation

1. Prepare qPCR Reactions MIX according to Table 4, and a 96-well plate layout template is shown in Table 5.

Table 4. qPCR Reaction MIX preparation

Tubes	Standard curver	NTC	NCS	Test sample
qPCR MIX	20 µL	20 µL	20 µL	20 µL
Samples	10 µL ST1 - ST5	10 µL DDB	10 µL purified NCS	10 µL purified test sample
Total Volume	30 µL	30 µL	30 µL	30 µL

Table 5. Example of 96-well plate layout

ST5	ST5	ST5							NTC	NTC	NTC	A
ST4	ST4	ST4							NCS	NCS	NCS	B
ST3	ST3	ST3										C
ST2	ST2	ST2										D
ST1	ST1	ST1										E
									S1	S1	S1	F
									S2	S2	S2	G
									S3	S3	S3	H
1	2	3	4	5	6	7	8	9	10	11	12	

- This example represents the assay for a standard curve with 5 concentration gradients (ST1 to ST5), 1 NTC, 1 NCS, 3 test samples (S1 to S3), and 3 replicates for each sample.

This is a watermark for the trial version, register to get the full one!

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.

Remove Watermark Now

■ qPCR program setting

Note: The following instructions apply only to the ABI7500 instrument with SDS v1.4. If you use a different instrument or software, refer to the applicable instrument or software documentation.

1. Create a new document, then in the Assay drop-down list, select **Standard Curve (Absolute Quantitation)**.
2. Click **New Detector**, then enter CAR or TCR in the Name field, select **FAM** in the Reporter Dye drop-down list and select **(none)** in the Quencher; Click **New Detector**, then enter SCG in the Name field, select **CY5** in the Reporter Dye drop-down list and select **(none)** in the Quencher; Click **New Detector**, then

enter IPC in the Name field. Select **VIC** in the Reporter Dye drop-down list and select **(none)** in the Quencher Dye drop-down list, then click **OK**.

3. Select **ROX** as the passive reference dye, then Click **Next**.
4. Select the applicable set of wells for the samples, then select the corresponding detector for each well.
5. Select Finish, and then set thermal-cycling conditions:
 - a. Set the thermal cycling reaction volume to 30 μL .
 - b. Set the temperature and time as following (Table 6):

Table 6. qPCR running program

Step	Temp.	Time(mm:sec)	Cycles
Activation	95°C	10:00	1
Denaturation	95°C	00:15	40
Annealing/extension	60°C*	01:00	

This is a watermark for the trial version, register to get the full one!

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.


Remove Watermark Now

Column drop-down list:

- a. NTC: target DNA detector task = NTC
 - b. NCS, test samples= Unknown
2. Set up the standard curve as shown in the following table 7:

Table 7. Settings for Standard curve

Tube label	Task	Conc. (copies/ μL)	
		Channel-FAM	Channel-CY5
ST1	Standard	2.33×10^6	2.33×10^6
ST2	Standard	2.33×10^5	2.33×10^5
ST3	Standard	2.33×10^4	2.33×10^4
ST4	Standard	2.33×10^3	2.33×10^3
ST5	Standard	2.33×10^2	2.33×10^2

3. Select the **Results** tab, then select Amplification Plot.
4. In the Data drop-down list, select **Delta Rn vs Cycle**.
5. In the Analysis Settings window, enter the following settings:
 - a. Select **Manual Ct**.
 - b. In the Threshold field, enter 0.02.
 - c. Select **Automatic Baseline**.
6. Click the button  in the toolbar, then wait the plate analyzing.
7. Select the **Result** tab > **Standard curve** tab, then verify the Slope, Intercept and R² values.
8. Select the Report tab, then achieve the mean quantity and standard deviation for each sample.
9. Select **File >> Export >> Results**. In the Save as type drop-down list, select **Results Export Files**, then click **Save**.

This is a watermark for the trial version, register to get the full one!

Benefits for registered users:

1. No watermark on the output documents.
2. Can operate scanned PDF files via OCR.
3. No page quantity limitations for converted PDF files.

Remove Watermark Now

12. Analyze the Ct value of IPC. Normally, the mean Ct-IPC value of the sample should be within ± 1.0 of the NCS Ct-IPC value. If the mean Ct-IPC value of the sample is significantly higher than the Ct-IPC value of the NCS, this indicates that the sample may be inhibitory to the assay.

13. Calculation of results:

CAR or TCR copies/cell = $2 \times$ CAR or TCR/SCG

Note: The parameter settings of the result analysis should be based on the specific model and the software version, and generally can also be automatically interpreted by the instrument.

Effective date: 08 Jul. 2024

This is a watermark for the trial version, register to get the full one!

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.

Remove Watermark Now

Support & Contact

SHENTEK

Huzhou Shenke Biotechnology Co., Ltd.

www.shentekbio.com

Address: 8th Floor, 6B Building, No.1366 Hongfeng Road, Huzhou313000, Zhejiang Province, China

E-mail: info@shentekbio.com

Phone: +1 (908) 822-3199 / (+86) 400-878-2189