

Anti-*E.coli* (K-12 & Alkaline Lysis) HCP Antibody

Product No.: 1301302-Ab01

■ Product information

Size	1 mg
Species Reactivity	<i>Escherichia coli</i> (<i>E. coli</i>)
Host/Isotype	Sheep / IgG
Class	Polyclonal
Type	Antibody
Conjugate	Unconjugated
Immunogen	<i>E.coli</i> HCPs (K-12 strains) produced chemically by alkaline lysis
Form	Liquid
Purification	Affinity Purified

Storage buffer

PBS

Contains

No preservatives

This is a watermark for the trial version, register to get the full one!

Storage

Maintain at 2-8°C for short-term use. For long-term storage, please store

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.

[Remove Watermark Now](#)

This product is for research use only and is not intended for diagnostic purposes.

■ Related product

Product Name	Product No.	Quantity
<i>E.coli</i> (K-12 & Alkaline Lysis) HCP ELISA	1301302	96 Tests

Effective date: 10 Jul. 2024

Support & Contact

SHENTEK

Huzhou Shenke Biotechnology Co., Ltd.

www.shentekbio.com

Address: 8th Floor, 6B Building, No.1366 Hongfeng Road, Huzhou313000, Zhejiang Province, China

E-mail: info@shentekbio.com

Phone: +1 (908) 822-3199 / (+86) 400-878-2189