

# Anti-*E.coli* (Protein Expression Strains) HCP Antibody

Product No.: 1301301-Ab01

## ■ Product information

<b>Size</b>	1 mg
<b>Species Reactivity</b>	<i>Escherichia coli</i> ( <i>E. coli</i> )
<b>Host/Isotype</b>	Sheep / IgG
<b>Class</b>	Polyclonal
<b>Type</b>	Antibody
<b>Conjugate</b>	Unconjugate
<b>Immunogen</b>	<i>E.coli</i> HCPs (BL21) produced by physical disruption
<b>Form</b>	Liquid
<b>Purification</b>	Affinity Purified

Storage buffer

PBS

Contains

No preservatives

This is a watermark for the trial version, register to get the full one!

Storage

Maintain at 2-8°C for short-term use. For long-term storage, please store

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.

[Remove Watermark Now](#)

## ■ Related product

Product Name	Product No.	Quantity
<i>E.coli</i> (Protein Expression Strains) HCP ELISA Kit (One-step ELISA)	1301301-1	96 Tests

Effective date: 10 Jul. 2024

## Support & Contact

**SHENTEK**

Huzhou Shenke Biotechnology Co., Ltd.

[www.shentekbio.com](http://www.shentekbio.com)

Address: 8th Floor, 6B Building, No.1366 Hongfeng Road, Huzhou313000, Zhejiang Province, China

E-mail: [info@shentekbio.com](mailto:info@shentekbio.com)

Phone: +1 (908) 822-3199 / (+86) 400-878-2189