

Anti-MDCK HCP Antibody

Product No.: 1301308-Ab01

■ Product Details

| | |
|---------------------------|---|
| Size | 1 mg |
| Species Reactivity | <i>Canis familiaris</i> , dog |
| Host/Isotype | Sheep / IgG |
| Class | Polyclonal |
| Type | Antibody |
| Conjugate | Unconjugated |
| Immunogen | MDCK HCPs produced by sonication cell lysis |
| Form | Liquid |
| Purification | Affinity Purified |

Buffer: PBS

Contains: No preservatives

This is a watermark for the trial version, register to get the full one!

Storage: Maintain at 2-8°C for short-term use. For long-term storage, please store

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.

[Remove Watermark Now](#)

This product is for research use only but not intended for diagnostic purposes.

■ Related Product

| Product Name | Product No. | Quantity |
|-------------------------------------|-------------|----------|
| MDCK HCP ELISA Kit (One-step ELISA) | 1301308 | 96 Tests |

Effective date: 10 Jul. 2024

Support & Contact

SHENTEK

Huzhou Shenke Biotechnology Co., Ltd.

www.shentekbio.com

Address: 8th Floor, 6B Building, No.1366 Hongfeng Road, Huzhou313000, Zhejiang Province, China

E-mail: info@shentekbio.com

Phone: +1 (908) 822-3199 / (+86) 400-878-2189