

Anti-Vero HCP Antibody

Product No.: 1301309-Ab01

■ Product information

Size	1 mg
Species Reactivity	African green monkey
Host/Isotype	Sheep / IgG
Class	Polyclonal
Type	Antibody
Conjugate	Unconjugated
Immunogen	Vero HCPs produced by sonication cell lysis
Form	Liquid
Purification	Affinity Purified

Storage buffer: PBS

Contains: No preservatives

This is a watermark for the trial version, register to get the full one!

Storage conditions: Maintain at 2-8°C for short-term use. For long-term storage, please store

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.

[Remove Watermark Now](#)

This product is for research use only and is not intended for diagnostic purposes.

■ Related product

Product Name	Product No.	Quantity
Vero HCP ELISA Kit (One-step ELISA)	1301309	96 Tests

Effective date: 10 Jul. 2024

Support & Contact

SHENTEK

Huzhou Shenke Biotechnology Co., Ltd.

www.shentekbio.com

Address: 8th Floor, 6B Building, No.1366 Hongfeng Road, Huzhou313000, Zhejiang Province, China

E-mail: info@shentekbio.com

Phone: +1 (908) 822-3199 / (+86) 400-878-2189