

Anti-CHO HCP Antibody

Product No.: 1301304-Ab01

■ Product information

Size	1 mg
Species Reactivity	Chinese hamster
Host/Isotype	Sheep / IgG
Class	Polyclonal
Type	Antibody
Conjugate	Unconjugated
Immunogen	CHO HCPs produced by CHO K1&CHO S Fed-batch supernatant
Form	Liquid
Purification	Affinity purified
Buffer	PBS
Contains	No preservatives
Applications	HCP coverage evaluations, Western Blot, ELISA
Storage	Maintain at 2-8°C for short-term use. For long-term storage, please store below -65°C.
Validity period	12 months

■ Attentions

This product is for research use only but not intended for diagnostic purposes.

■ Related product

Product Name	Product No.	Quantity
CHO HCP ELISA Kit (One-step ELISA)	1301304-1	96 Tests

Effective date: 10 Jul. 2024

Support & Contact



Huzhou Shenke Biotechnology Co., Ltd.

www.shentekbio.com

Address: 8th Floor, 6B Building, No.1366 Hongfeng Road, Huzhou313000, Zhejiang Province, China

E-mail: info@shentekbio.com

Phone: +1 (908) 822-3199 / (+86) 400-878-2189