

Anti-*E.coli* (Protein Expression Strains) HCP Antibody

Product No.: 1301301-Ab01

■ Product information

Size	1 mg
Species Reactivity	<i>Escherichia coli</i> (<i>E. coli</i>)
Host/Isotype	Sheep / IgG
Class	Polyclonal
Type	Antibody
Conjugate	Unconjugate
Immunogen	<i>E.coli</i> HCPs (BL21) produced by physical disruption
Form	Liquid
Purification	Affinity Purified
Storage buffer	PBS
Contains	No preservatives
Applications	HCP Coverage evaluations, Western Blot, ELISA
Storage	Maintain at 2-8°C for short-term use. For long-term storage, please store below -65°C.
Validity period	12 months

■ Attentions

This product is for research use only and is not intended for diagnostic purposes.

■ Related product

Product Name	Product No.	Quantity
<i>E.coli</i> (Protein Expression Strains) HCP ELISA Kit (One-step ELISA)	1301301-1	96 Tests

Effective date: 10 Jul. 2024

Support & Contact



Huzhou Shenke Biotechnology Co., Ltd.

www.shentekbio.com

Address: 8th Floor, 6B Building, No.1366 Hongfeng Road, Huzhou313000, Zhejiang Province, China

E-mail: info@shentekbio.com

Phone: +1 (908) 822-3199 / (+86) 400-878-2189