

Anti-HEK293 HCP Antibody

Product No.: 1301311-Ab01

■ Product information

Size	1 mg
Species Reactivity	Human
Host/Isotype	Sheep / IgG
Class	Polyclonal
Type	Antibody
Conjugate	Unconjugated
Immunogen	HEK293 HCPs produced by sonication cell lysis
Form	Liquid
Purification	Affinity Purified
Storage buffer	PBS
Contains	No preservatives
Applications	HCP Coverage evaluations, Western Blot, ELISA
Storage conditions	Maintain at 2-8°C for short-term use. For long-term storage, please store below -65°C.
Validity period	12 months

■ Attentions

This product is for research use only and is not intended for diagnostic purposes.

■ Related product

Product Name	Product No.	Quantity
HEK293 HCP ELISA Kit (One-step ELISA)	1301311	96 Tests

Effective date: 10 Jul. 2024

Support & Contact

SHENTEK

Huzhou Shenke Biotechnology Co., Ltd.

www.shentekbio.com

Address: 8th Floor, 6B Building, No.1366 Hongfeng Road, Huzhou313000, Zhejiang Province, China

E-mail: info@shentekbio.com

Phone: +1 (908) 822-3199 / (+86) 400-878-2189